# Culture Medium for Activation and Expansion of Activated Human NK Cell

# KBM 501 & 502/502B

#### These products are registered in MF and the FDA Drug Master File.

#### **Overview**


KBM501 is optimized for activation of human peripheral blood natural killer cells. KBM502 and KBM502B are used for the subsequent expansion culture of the activated lymphocytes with KBM501.

### Characteristics

- No proteins other than human serum albumin, human transferase, recombinant human insulin and recombinant human interleukin 2 have been added.
- Capable of minimize pH variation by enhanced buffer capacity.
- Kanamycin sulphate is contained as an antibiotic.
- The bag of KBM502B is made of a material with excellent gas permeability and its safety and quality have been confirmed based on the plastic drug container test method of the Japanese Pharmacopoeia.

## Cell Culture Example

Proliferative response of IL2-dependent human lymphocytes on immobilized anti-CD3 antibody



#### [Culture Protocols]

Seeded Cells :	Normal human peripheral blood monocytes				
Medium :	KBM 501(recombinant human interleukin 2 : 2810 IU/mL)				
	1% inactivated autoplasma was added during the culture.				
Culture Vessel : 12-well plate on which anti-human-CD3 antibody is immobilized.					
Protocols :	Cultured at a seeding density of $2.7 \times 10^5$ cells/mL, under static condition at $37^{\circ}$ C • 5% CO2. Expand the scale of culture depending on cell proliferation with medium containing plasma and recombinant human interleukin-2.				
Analysis :	At day 10, surface markers of activated lymphocytes were				

analyzed by flow cytometry, and cytotoxicity to K562 cells were measured with TeraScan.

**KBM 502** 

ドトアルフ 細胞カナマ


ロット美号:K2M108 本計画目:2011/03

3-22

Cytotoxicity of activated lymphocytes to K562 cells

**KBM 501** 

MRA+1



Surface markers of activated lymphocytes

		Incidence (%)			
	Surface marker	$\begin{array}{c} \text{KBM} \\ 501 \end{array}$	Medium A	Medium A5	
-	CD3- CD56+	1.74	1.23	0.54	
	CD3+ CD56+	15.52	32.93	35.91	
	CD3- CD56-	0.93	2.06	0.68	
	CD3+ CD56-	81.81	63.78	62.87	

Product No.	Product Name	Size	Price	Shelf time	Storage
16025015	KBM 501 (IL-2: 2810 IU/mL)	500  mL	JPY 20,000	8 months	2-8°C
16025020	KBM 502 (IL-2: 281 IU/mL)	1000 mL	JPY 16,000	8 months	2-8°C
1602502B	KBM 502B (IL-2: 281 IU/mL)	1000 mL	JPY 17,500	8 months	2-8°C

\* This product is sold for research purpose only.