



COMPOSITION

Dulbecco's Phosphate Buffered Saline DPBS

Inorganic Salts

	CAS number	L0615 Liquid mg/l	L0616 Liquid mg/l	X0515 Liquid mg/l	X0520 Liquid mg/l	P0750 Liquid mg/l
Calcium Chloride Dihydrate	10035-04-8	/	/	/	1330	/
Magnesium Chloride Hexahydrate	7791-18-6	/	/	/	1000	/
Potassium Chloride	7447-40-7	200	/	2000	2000	200
Potassium Phosphate Monobasic Anhydrous	7778-77-0	200	144	2000	2000	200
Sodium Chloride	7647-14-5	8000	9000	80000	80000	8000
Sodium Phosphate Dibasic Anhydrous	7558-79-4	1150	795	11500	11500	1150