



# COMPOSITION

## Hank's Balanced Salts Solution

	CAS number	L0605 Liquid mg/l	L0606 Liquid mg/l	L0607 Liquid mg/l	L0611 Liquid mg/l	L0612 Liquid mg/l
Calcium Chloride Dihydrate	10035-04-8	/	185	/	/	185
Magnesium Sulfate Anhydrous	7487-88-9	/	97,67	/	/	97,67
Potassium Chloride	7447-40-7	400	400	400	400	400
Potassium Phosphate Monobasic Anhydrous	7778-77-0	60	60	60	60	60
Sodium Bicarbonate	144-55-8	/	350	350	350	350
Sodium Chloride	7647-14-5	8000	8000	8000	8000	8000
Sodium Phosphate Dibasic Anhydrous	7558-79-4	47,88	47,88	47,88	47,88	47,88
D-Glucose Anhydrous	50-99-7	1000	1000	1000	1000	1000
Phenol Red Sodium Salt	34487-61-1	/	11	/	11	/