



## Technical data sheet

Ref : FT.L0093an

Page : 1/2

Version date : 25/07/2022

### DMEM Ham's F12 w/ L-Glutamine w/ 15mM Hepes

**CAT N°:** L0093

**Storage conditions:** +2°C to +8°C in the dark

**Shelf life:** 12 months

**Composition:** Displayed on website and in catalog; also available on request.

**Colour:** Clear red solution

**pH:** 7.3 ± 0.3

**Osmolality:** 299 mOsm/kg ± 10 %

**Endotoxin:** < 1 EU/ml

**Sterility tests:**

- Bacteria in aerobic and anaerobic conditions
- Fungi and yeasts

**Cell Growth test:**

Medium tested for the ability to support Hela cell growth.

**Other tests:** Not applicable

**Recommended use:**

- Respect storage conditions of the product
- Do not use the product after its expiry date
- Store product in an area protected from light (not necessary for saline solutions).
- Manipulate the product in aseptic conditions (e.g.: under laminar air flow)
- Wear clothes adapted to the manipulation of the product to avoid contamination (e.g.: gloves, mask, hygiene cap, overall...)

The product is intended to be used in vitro for research or further manufacturing only and not for use as an Active Pharmaceutical Ingredient or food or animal feed.

**Application:**

Studies to determine the nutritional requirements of many cells have been in progress since Eagle's first reports. The major essential nutrients were identified and work became focused on the media requirements of individual cell types. Many media designed for these purposes are now available. Among the first of these media, developed initially to study hormonal requirements of cells in culture, was a mixture of DMEM medium and Ham's F12 medium, known as DMEM Ham's F12.

**Uses:**

Supplements, such as antibiotics, should be added as sterile supplements to the medium. Storage conditions and shelf-life of supplemented product will be affected by the nature of the supplements.

**Signs of Deterioration:**



## Technical data sheet

Ref : FT.L0093an  
Page : 2/2

Version date : 25/07/2022

Medium should be clear and free of particulate and flocculent material. Do not use if medium is cloudy or contains precipitate.

Other evidence of deterioration may include colour change or degradation of physical or performance characteristics.

**Remarks:** Not applicable