

COMPOSITION

MEM Vitamins

	CAS number	X0556 Liquid mg/l
Sodium Chloride	7647-14-5	8500
Choline Chloride	67-48-1	100
D-Ca Pantothenate	137-08-6	100
Folic Acid	59-30-3	100
Myo-Inositol	87-89-8	200
Nicotinamide (Nicotinic acid amide)	98-92-0	100
Pyridoxal Hydrochloride	65-22-5	100
Riboflavin	83-88-5	10
Thiamine Hydrochloride	67-03-8	100

MEM Non Essential Amino Acids

	CAS number	X0557 Liquid mg/l
Glycine	56-40-6	750
L-Alanine	56-41-7	890
L-Asparagine Monohydrate	5794-13-8	1500
L-Aspartic Acid	56-84-8	1330
L-Glutamic Acid	56-86-0	1470
L-Proline	147-85-3	1150
L-Serine	56-45-1	1050