




FETAL BOVINE SERUM IgG Depleted

YOUR CELLS WILL FEEL THE DIFFERENCE

Some experiments require low levels of specific antibody. Fetal bovine serum (FBS) has natural levels of IgG that are too high for some cell culture applications. Biowest can offer you a perfect alternative without changing the performance of FBS.

- Performance is equivalent to non-treated FBS
- Ultra low IgG levels < 5 µg/ml
- No modification on the cell growth promotion
- No modification on the cell viability
- Ideal solution for cell culture and protein purification applications where naturally occurring levels of IgG are too high



FBS IgG depleted has been treated with proprietary chromatography method. This treatment reduces the level of IgG.

Ordering Information

Cat N°	Product Name
S1811-100	FBS South America, IgG depleted, 100mL
S1811-500	FBS South America, IgG depleted, 500mL

Please visit www.biowest.net for more information

www.biowest.net

2, rue du Vieux Bourg, 49340 NUAILLÉ
Tel.: (+33) 241 464 242
Email: biowest@biowest.net